

ALK

Catalog Number: 21280

Gene Symbol: ALK, CD246, NBLST3

Description: Anti-ALK Rabbit Polyclonal Antibody

Background: Anaplastic lymphoma kinase (ALK) also known as ALK tyrosine kinase receptor or CD246 (cluster of differentiation 246) is a neuronal orphan receptor tyrosine kinase. ALK is essentially and transiently expressed in specific regions of the central and peripheral nervous systems and plays an important role in the genesis and differentiation of the nervous system. Abnormal ALK, by fusion with other genes, gaining additional gene copies or mutations, can be oncogenic.

Immunogen: Purified ALK protein, human origin.

Tested applications: ELISA, WB, IHC

Recommended dilutions:

ELISA: 1:1000-1:5000

WB: 1:100-1:1000

IHC: 1:50-1:200

Concentration: 0.2mg/ml

Host:Rabbit

Clonality: Polyclonal

Purity: Purified from serum

Format: Liquid

Storage buffer:


Preservative: no

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 50% glycerol

Species Reactivity: Recognizes ALK of vertebrates.

Storage Conditions: Store at -20°C. Avoid freeze / thaw cycles

Western blot:



Western blot analysis of recombinant ALK protein. Purified His-tagged ALK protein was blotted with anti-ALK rabbit polyclonal antibody (Cat. #21280).

For research use only. Not for diagnostic or therapeutic applications.


24 Whitewoods Lane, Malvern, PA 19355

E-mail: sale@neweastbio.com

Tel: (610) 945-2007 FAX: (610) 945-2008

Web: www.neweastbio.com

Immunohistochemistry:



Immunohistochemical analysis of paraffin-embedded human thyroid cancer tissue with anti-ALK polyclonal antibody (Cat. #21280). Tissue samples were fixed with formaldehyde and blocked with 1% serum for 15 min at 37 °C. Antigen retrieval was by heat mediation in citrate buffer (pH6). Samples were then incubated with primary antibody (1:50 dilution) overnight at 4°C. A HRP-conjugated Goat anti- Rabbit IgG (1:50 dilution) was used as secondary antibody.