

Catalog #	Product name	Size / Price (s)
26261	Anti-RAS (A146T) Mouse Monoclonal Antibody	100ul / \$495
26270	Anti-MEK1 (N382H) Mouse Monoclonal Antibody	100ul / \$495
26305	Anti-CTNNB1 (T41A) Mouse Monoclonal Antibody	100ul / \$495
26242	Anti-EML4-ALK Mouse Monoclonal Antibody	100ul / \$495
26168	Anti-CTNNB1 (S37F) Mouse Monoclonal Antibody	100ul / \$495
26216	Anti-PIK3CA (H1047Y) Mouse Monoclonal Antibody	100ul / \$495
26210	Anti-PIK3CA (E542Q) Mouse Monoclonal Antibody	100ul / \$495
26215	Anti-PIK3CA (D549N) Mouse Monoclonal Antibody	100ul / \$495
26263	Anti-MEK1 (F53L) Mouse Monoclonal Antibody	100ul / \$495
26266	Anti-MEK1 (P264S) Mouse Monoclonal Antibody	100ul / \$495
26303	Anti-CTNNB1 (D32N) Mouse Monoclonal Antibody	100ul / \$495
26304	Anti-CTNNB1 (S33F) Mouse Monoclonal Antibody	100ul / \$495
26306	Anti-CTNNB1 (T41I) Mouse Monoclonal Antibody	100ul / \$495
26296	Anti-CTNNB1 (S37Y) Mouse Monoclonal Antibody	100ul / \$495
26251	Anti-FGFR1 (R78H) Mouse Monoclonal Antibody	100ul / \$495
26253	Anti-FGFR1 (P252T) Mouse Monoclonal Antibody	100ul / \$495
26300	Anti-CTNNB1 (G34E) Mouse Monoclonal Antibody	100ul / \$495
26308	Anti-CTNNB1 (S45F) Mouse Monoclonal Antibody	100ul / \$495
26294	Anti-FGFR1 Mouse Monoclonal Antibody	100ul / \$495
26259	Anti-FGFR1 (V664L) Mouse Monoclonal Antibody	100ul / \$495
26252	Anti-FGFR1 (S125L) Mouse Monoclonal Antibody	100ul / \$495
26257	Anti-FGFR1 (R576W) Mouse Monoclonal Antibody	100ul / \$495

26302	Anti-CTNNB1 (D32Y) Mouse Monoclonal Antibody	100ul / \$495
26309	Anti-CTNNB1 (S45Y) Mouse Monoclonal Antibody	100ul / \$495
26287	Anti-EGFR (L747-P753) S Mouse Monoclonal Antibody	100ul / \$495
26254	Anti-FGFR1 (P252S) Mouse Monoclonal Antibody	100ul / \$495
26299	Anti-CTNNB1 (G34V) Mouse Monoclonal Antibody	100ul / \$495
26314	Anti-JAK2 (R683G) Mouse Monoclonal Antibody	100ul / \$495
26315	Anti-JAK2 (R683S) Mouse Monoclonal Antibody	100ul / \$495
26319	Anti-BRCA1 (C61G) Mouse Monoclonal Antibody	100ul / \$495
26328	Anti-GNAQ (Q209L) Mouse Monoclonal Antibody	100ul / \$495
26330	Anti-GNAQ (Q209R) Mouse Monoclonal Antibody	100ul / \$495
26332	Anti-NFE2L2 (L30F) Mouse Monoclonal Antibody	100ul / \$495
26336	Anti-NFE2L2 (E79Q) Mouse Monoclonal Antibody	100ul / \$495
26338	Anti-RHEB (Y35N) Mouse Monoclonal Antibody	100ul / \$495
26167	Anti- β -catenin (E53K) Mouse Monoclonal Antibody	100ul / \$495
26169	Anti- β -catenin (G38D) Mouse Monoclonal Antibody	100ul / \$495
26166	Anti- β -catenin (G48D) Mouse Monoclonal Antibody	100ul / \$495
26170	Anti- β -catenin (A43V) Mouse Monoclonal Antibody	100ul / \$495
26154	Anti-Abl (Y253H) Mouse Monoclonal Antibody	100ul / \$495
26137	Anti-AKT1 (E49K) Mouse Monoclonal Antibody	100ul / \$495
26226	Anti-ALK (C1156Y) Mouse Monoclonal Antibody	100ul / \$495
26151	Anti-ALK (R1275Q) Mouse Monoclonal Antibody	100ul / \$495
26141	Anti-ATM (F858L) Mouse Monoclonal Antibody	100ul / \$495
26142	Anti-ATM (Q2442P) Mouse Monoclonal Antibody	100ul / \$495

26143	Anti-ATM (R2443Q) Mouse Monoclonal Antibody	100ul / \$495
26140	Anti-ATM (R3008H) Mouse Monoclonal Antibody	100ul / \$495
26109	Anti-B-Raf (N581S) Mouse Monoclonal Antibody	100ul / \$495
26184	Anti-B-Raf (V600D) Mouse Monoclonal Antibody	100ul / \$495
26039	Anti-B-Raf(V600E) Mouse Monoclonal Antibody	100ul / \$495
26110	Anti-B-Raf (V600K) Mouse Monoclonal Antibody	100ul / \$495
26182	Anti-B-Raf (V600R) Mouse Monoclonal Antibody	100ul / \$495
26196	Anti-BRAF (G469E) Mouse Monoclonal Antibody	100ul / \$495
26117	Anti-EGFR (A839T) Mouse Monoclonal Antibody	100ul / \$495
26120	Anti-EGFR (D761Y) Mouse Monoclonal Antibody	100ul / \$495
26115	Anti-EGFR(G719A) Mouse Monoclonal Antibody	100ul / \$495
26083	Anti-EGFR(G719S) Mouse Monoclonal Antibody	100ul / \$495
26119	Anti-EGFR (G863D) Mouse Monoclonal Antibody	100ul / \$495
26221	Anti-EGFR (L858Q) Mouse Monoclonal Antibody	100ul / \$495
26082	Anti-EGFR (L858R) Mouse Monoclonal Antibody	100ul / \$495
26077	Anti-EGFR (L861Q) Mouse Monoclonal Antibody	100ul / \$495
26116	Anti-EGFR (N826S) Mouse Monoclonal Antibody	100ul / \$495
26121	Anti-EGFR (V769L) Mouse Monoclonal Antibody	100ul / \$495
26176	Anti-EGFR (vIII) Mouse Monoclonal Antibody	100ul / \$495
26080	Anti-ErbB2 (V777A) mouse monoclonal Antibody	100ul / \$495
26249	Anti-FGFR1 (P252R) Mouse Monoclonal Antibody	100ul / \$495
26386	Anti-FGFR1 N-Terminal Mouse Monoclonal Antibody	100ul / \$495
26329	Anti-GNAQ (Q209P) Mouse Monoclonal Antibody	100ul / \$495

26081	Anti-IDH1 (R132H) Mouse Monoclonal Antibody	100ul / \$495
26406	Anti-IDH1 (R132P) Mouse Monoclonal Antibody	100ul / \$495
26160	Anti-IDH1 (R132S) Mouse Monoclonal Antibody	100ul / \$495
26165	Anti-IDH2 (R140Q) Mouse Monoclonal Antibody	100ul / \$495
26084	Anti-IDH2 (R172H) Mouse Monoclonal Antibody	100ul / \$495
26163	Anti-IDH2 (R172K) Mouse Monoclonal Antibody	100ul / \$495
26408	Anti-IDH2 (R172S) Mouse Monoclonal Antibody	100ul / \$495
26164	Anti-IDH2 (R172W) Mouse Monoclonal Antibody	100ul / \$495
26264	Anti-MEK1 (I111S) Mouse Monoclonal Antibody	100ul / \$495
26222	Anti-MEK1 (Q56P) Mouse Monoclonal Antibody	100ul / \$495
26331	Anti-NFE2L2 (D29G) Mouse Monoclonal Antibody	100ul / \$495
26337	Anti-NFE2L2 (E79K) Mouse Monoclonal Antibody	100ul / \$495
21145	Anti-PIK3CA(E545K) Rabbit Polyclonal Antibody	100ul / \$495
26216	Anti-PIK3CA (H1047Y) Mouse Monoclonal Antibody	100ul / \$495
26138	Anti-PIK3CA (Q546K) Mouse Monoclonal Antibody	100ul / \$495
21146	Anti-PIK3CA(H1047R) Rabbit Polyclonal Antibody	100ul / \$495
26186	Anti-RAS (G12S) Mouse Monoclonal Antibody	100ul / \$495
26191	Anti-HRAS (G13A) Mouse Monoclonal Antibody	100ul / \$495
26195	Anti-Ras (Q61H) Mouse Monoclonal Antibody	100ul / \$495
26193	Anti-Ras (Q61L) Mouse Monoclonal Antibody	100ul / \$495
26192	Anti-Ras (Q61R) Mouse Monoclonal Antibody	100ul / \$495
26036	Anti-Ras(G12D) Mouse Monoclonal Antibody	100ul / \$495
26038	Anti-Ras (G13D) Mouse Monoclonal Antibody	100ul / \$495

26422	Anti-SMAD4 (R361H) Mouse Monoclonal Antibody	100ul / \$495
26158	Anti-STK11 Mouse Monoclonal Antibody	100ul / \$495
26072	Anti-TP53(R175H) Mouse Monoclonal Antibody	100ul / \$495
21165	Anti-TP53(R282W) Rabbit Polyclonal Antibody	100ul / \$495
21161	Anti-TP53(R175H) Rabbit Polyclonal Antibody	100ul / \$495
26076	Anti-TP53(R282W) Mouse Monoclonal Antibody	100ul / \$495
26290	Anti-ALK (D1225N) Mouse Monoclonal Antibody	100ul / \$495
26200	Anti-BRAF (D594G) Mouse Monoclonal Antibody	100ul / \$495
26347	Anti-FBXW7 (R479L) Mouse Monoclonal Antibody	100ul / \$495
26268	Anti-MEK1 (E203K) Mouse Monoclonal Antibody	100ul / \$495
26281	Anti-KIT (V560D) Mouse Monoclonal Antibody	100ul / \$495
26272	Anti-KIT (L576P) Mouse Monoclonal Antibody	100ul / \$495
26273	Anti-KIT (H697Y) Mouse Monoclonal Antibody	100ul / \$495
26275	Anti-KIT (W557R) Mouse Monoclonal Antibody	100ul / \$495
26277	Anti-KIT (V559D) Mouse Monoclonal Antibody	100ul / \$495
26278	Anti-KIT (K642E) Mouse Monoclonal Antibody	100ul / \$495
26219	Anti-PIK3CA (G1049S) Mouse Monoclonal Antibody	100ul / \$495
26234	Anti-PIK3CA (E542K) Mouse Monoclonal Antibody	100ul / \$495
26214	Anti-PIK3CA (Q546R) Mouse Monoclonal Antibody	100ul / \$495
26213	Anti-PIK3CA (Q546P) Mouse Monoclonal Antibody	100ul / \$495
26235	Anti-EGFR (E746-A750)del Mouse Monoclonal Antibody	100ul / \$495
26199	Anti-BRAF (V600A) Mouse Monoclonal Antibody	100ul / \$495
26462	Anti-BRAF (V600M) Mouse Monoclonal Antibody	100ul / \$495

26351-1	Anti-BRCA1 (1-120) Mouse Monoclonal Antibody	100ul / \$495
26351-2	Anti-BRCA1 (1720-1863) Mouse Monoclonal Antibody	100ul / \$495
26364-2	Anti-BRCA2 (1475-1748) Mouse Monoclonal Antibody	100ul / \$495
26364-3	Anti-BRCA2 (2213-2543) Mouse Monoclonal Antibody	100ul / \$495
26364-4	Anti-BRCA2 (2793-3030) Mouse Monoclonal Antibody	100ul / \$495
26364-1	Anti-BRCA2 (352-643) Mouse Monoclonal Antibody	100ul / \$495
26324	Anti-CTCF (K344E) Mouse Monoclonal Antibody	100ul / \$495
26326	Anti-CTCF (R339W) Mouse Monoclonal Antibody	100ul / \$495
26307	Anti-CTNNB1 (S45P) Mouse Monoclonal Antibody	100ul / \$495
26483	Anti-EGFR (S768I) Mouse Monoclonal Antibody	100ul / \$495
26481	Anti-EGFR (T854A) Mouse Monoclonal Antibody	100ul / \$495
26382	Anti-ERBB2 (D769Y) Mouse Monoclonal Antibody	100ul / \$495
26231	Anti-IDH2 (R172G) Mouse Monoclonal Antibody	100ul / \$495
26274	Anti-KIT (D820E) Mouse Monoclonal Antibody	100ul / \$495
26288	Anti-KIT (P577-D579 DEL) Mouse Monoclonal Antibody	100ul / \$495
26289	Anti-KIT (V560DEL) Mouse Monoclonal Antibody	100ul / \$495
26335	Anti-NFE2L2 (D77V) Mouse Monoclonal Antibody	100ul / \$495
26208	Anti-PIK3CA (E545G) Mouse Monoclonal Antibody	100ul / \$495
26216	Anti-PIK3CA (H1047Y) Mouse Monoclonal Antibody	100ul / \$495
26478	Anti-RAS (A146P) Mouse Monoclonal Antibody	100ul / \$495
26479	Anti-RAS (A146V) Mouse Monoclonal Antibody	100ul / \$495
26474	Anti-RAS (G13R) Mouse Monoclonal Antibody	100ul / \$495
26477	Anti-RAS (G13V) Mouse Monoclonal Antibody	100ul / \$495

26468	Anti-RAS (K117N) Mouse Monoclonal Antibody	100ul / \$495
26146	Anti-RET (C634Y) Mouse Monoclonal Antibody	100ul / \$495
26425	Anti-SMO (A68V) Mouse Monoclonal Antibody	100ul / \$495
26444	Anti-NOTCH1 (M2011R) Mouse Monoclonal Antibody	100ul / \$495
26901	Anti-Active G α i-GTP Mouse Monoclonal Antibody	30ul / \$495
26902	Anti-Active G α 13-GTP Mouse Monoclonal Antibody	30ul / \$495
26903	Anti-Active Rac1-GTP Mouse Monoclonal Antibody	30ul / \$495
26904	Anti-Active RhoA-GTP Mouse Monoclonal Antibody	30ul / \$495
26905	Anti-Active Cdc42-GTP Mouse Monoclonal Antibody	30ul / \$495
26906	Anti-Active G α s-GTP Mouse Monoclonal Antibody	30ul / \$495
26907	Anti-Active G α o-GTP Mouse Monoclonal Antibody	30ul / \$495
26908	Anti-Active G α z-GTP Mouse Monoclonal Antibody	30ul / \$495
26909	Anti-Active Ras-GTP Mouse Monoclonal Antibody	30ul / \$495
26910	Anti-Active Rheb-GTP Mouse Monoclonal Antibody	30ul / \$495
26911	Anti-Active Rab5-GTP Mouse Monoclonal Antibody	30ul / \$495
26912	Anti-Active Rap1-GTP Mouse Monoclonal Antibody	30ul / \$495
26913	Anti-Active Ral-GTP Mouse Monoclonal Antibody	30ul / \$495
26915	Anti-Active Ran-GTP Mouse Monoclonal Antibody	30ul / \$495
26916	Anti-Active Sar1-GTP Mouse Monoclonal Antibody	30ul / \$495
26918	Anti-Active Arf6-GTP Mouse Monoclonal Antibody	30ul / \$495
26919	Anti-Active Rab11-GTP Mouse Monoclonal Antibody	30ul / \$495
26920	Anti-Active Rab3-GTP Mouse Monoclonal Antibody	30ul / \$495
26921	Anti-Active Rab4-GTP Mouse Monoclonal Antibody	30ul / \$495

26922	Anti-Active Rab35-GTP Mouse Monoclonal Antibody	30ul / \$495
26923	Anti-Active Rab7-GTP Mouse Monoclonal Antibody	30ul / \$495
26924	Anti-Active Arl1-GTP Mouse Monoclonal Antibody	30ul / \$495
26925	Anti-Active Arl3-GTP Mouse Monoclonal Antibody	30ul / \$495
21001	Anti-Orai1-L1 Rabbit Polyclonal Antibody	100ul / \$219
21002	Anti-Orai1-L2 Rabbit Polyclonal Antibody	100ul / \$219
21011	Anti- β -Actin Rabbit Polyclonal Antibody	100ul / \$219
21012	Anti-GAPDH Rabbit Polyclonal Antibody	100ul / \$219
21089	Anti- β -Tubulin Rabbit Polyclonal Antibody	100ul / \$219
21108	Anti-Cystatin C Rabbit Polyclonal Antibody	100ul / \$219
21110	Anti-Kim-1 Rabbit Polyclonal Antibody	100ul / \$219
21112	Anti-Albumin Rabbit Polyclonal Antibody	100ul / \$219
21116	Anti-Transferrin Rabbit Polyclonal Antibody	100ul / \$219
21118	Anti-SLC7A5 Rabbit Polyclonal Antibody	100ul / \$219
21123	Anti-MUC1 Rabbit Polyclonal Antibody	100ul / \$219
21126	Anti-Septin9 Rabbit Polyclonal Antibody	100ul / \$219
21128	Anti-EYA2 Rabbit Polyclonal Antibody	100ul / \$219
21236	Anti-SUMO-1 Rabbit Polyclonal Antibody	100ul / \$219
21237	Anti-SUMO-2/3 Rabbit Polyclonal Antibody	100ul / \$219
26012	Anti-Fascin(N-terminal peptide) Mouse Monoclonal Antibody	100ul / \$239
26013	Anti-Fascin(internal peptide) Mouse Monoclonal Antibody	100ul / \$239
26023	Anti-Apolipoprotein A1 Mouse Monoclonal Antibody	100ul / \$239
26025	Anti-Apolipoprotein A1 Mouse Monoclonal Antibody	100ul / \$239

26043	Anti-Clusterin Mouse Monoclonal Antibody	100ul / \$239
26044	Anti-Cystatin C Mouse Monoclonal Antibody	100ul / \$239
26047	Anti-TFF3 Mouse Monoclonal Antibody	100ul / \$239
26050	Anti-Transferrin Mouse Monoclonal Antibody	100ul / \$239
26051	Anti-BMP3 Mouse Monoclonal Antibody	100ul / \$239
26052	Anti-EYA2 Mouse Monoclonal Antibody	100ul / \$239
26054	Anti-Vimentin Mouse Monoclonal Antibody	100ul / \$239
26071	Anti-Ric-8A Mouse Monoclonal Antibody	100ul / \$239
26085	Anti-Cystatin C Mouse Monoclonal Antibody	100ul / \$239
26090	Anti-Albumin Mouse Monoclonal Antibody	100ul / \$239
26096	Anti-EYA2 Mouse Monoclonal Antibody	100ul / \$239
26131	Anti-Survivin Mouse Monoclonal Antibody	100ul / \$239
26133	Anti-GRB7 Mouse Monoclonal Antibody	100ul / \$239
26136	Anti-IDH1 Mouse Monoclonal Antibody	100ul / \$239
82001	Anti-Fascin ELISA Kit	100ul / \$495
83901	Anti-Human Apolipoprotein A1 EIA Kit	100ul / \$385
84601	Anti-Cystatin C ELISA Kit	100ul / \$385
26001-1	Anti-cGMP Mouse Monoclonal Antibody	100ul / \$330
26002-1	Anti-cAMP Mouse Monoclonal Antibody	100ul / \$330
26003	Anti-Gai Mouse Monoclonal Antibody	100ul / \$239
26004	Anti-Ga13 Mouse Monoclonal Antibody	100ul / \$239
26005	Anti-Rac Mouse Monoclonal Antibody	100ul / \$239
26006	Anti-Gas Mouse Monoclonal Antibody	100ul / \$239

26007	Anti-RhoA Mouse Monoclonal Antibody	100ul / \$239
26008	Anti-Cdc42 Mouse Monoclonal Antibody	100ul / \$239
26009	Anti-Gao Mouse Monoclonal Antibody	100ul / \$239
26011	Anti-Gaz Mouse Monoclonal Antibody	100ul / \$239
26014	Anti-Ras Mouse Monoclonal Antibody	100ul / \$239
26015	Anti-Rheb Mouse Monoclonal Antibody	100ul / \$239
26016	Anti-Rab5 Mouse Monoclonal Antibody	100ul / \$239
26018	Anti-RalB Mouse Monoclonal Antibody	100ul / \$239
26019	Anti-Arf1 Mouse Monoclonal Antibody	100ul / \$239
26020	Anti-Ran Mouse Monoclonal Antibody	100ul / \$239
26021	Anti-Sar1a Mouse Monoclonal Antibody	100ul / \$239
26023	Anti-Apolipoprotein A1 Mouse Monoclonal Antibody	100ul / \$239
26025	Anti-Apolipoprotein A1 Mouse Monoclonal Antibody	100ul / \$239
26026	Anti-Ga12 Mouse Monoclonal Antibody	100ul / \$239
26027	Anti-Ga14 Mouse Monoclonal Antibody	100ul / \$239
26029	Anti-Arf6 Mouse Monoclonal Antibody	100ul / \$239
26030	Anti-Rab11 Mouse Monoclonal Antibody	100ul / \$239
26031	Anti-Rab3A Mouse Monoclonal Antibody	100ul / \$239
26033	Anti-Rab35 Mouse Monoclonal Antibody	100ul / \$239
26055	Anti-Rab7 Mouse Monoclonal Antibody	100ul / \$239
26056	Anti-Arl1 Mouse Monoclonal Antibody	100ul / \$239
26060	Anti-Gq family Mouse Monoclonal Antibody	100ul / \$239
26064	Anti-Rac3 Mouse Monoclonal Antibody	100ul / \$239

26065	Anti-Ral Mouse Monoclonal Antibody	100ul / \$239
26070	Anti-Arl3 Mouse Monoclonal antibody	100ul / \$239
26135	Anti-Bag1 Mouse Monoclonal Antibody	100ul / \$239
26294	Anti-FGFR1 Mouse Monoclonal Antibody	100ul / \$239
26356	Anti-GNAQ Mouse Monoclonal Antibody	100ul / \$239
28004	Anti-His-tag Mouse Monoclonal Antibody	100ul / \$239
26358	Anti-NFE2L2 Mouse Monoclonal Antibody	100ul / \$239
26354	Anti-CTCF Mouse Monoclonal Antibody	100ul / \$239
26091	Anti-Vimentin Mouse Monoclonal Antibody	100ul / \$239
21003	Anti-Rac1 Rabbit Polyclonal Antibody	100ul / \$219
21005	Anti-Ga13 Rabbit Polyclonal Antibody	100ul / \$219
21007	Anti-Gas Rabbit Polyclonal Antibody	100ul / \$219
21010	Anti-Cdc42 Rabbit Polyclonal Antibody	100ul / \$219
21015	Anti-Gao Rabbit Polyclonal Antibody	100ul / \$219
21016	Anti-Gaz Rabbit Polyclonal Antibody	100ul / \$219
21017	Anti-RhoA Rabbit Polyclonal Antibody	100ul / \$219
21019	Anti-Cdc42 Chicken Polyclonal Antibody	100ul / \$219
21020	Anti-Ga13 Rabbit Polyclonal Antibody	100ul / \$219
21021	Anti-Ras Rabbit Polyclonal Antibody	100ul / \$219
21030	Anti-Arf1 Rabbit Polyclonal Antibody	100ul / \$219
21032	Anti-Arf6 Rabbit Polyclonal Antibody	100ul / \$219
21033	Anti-RalB Rabbit Polyclonal Antibody	100ul / \$219
21034	Anti-RalA Rabbit Polyclonal Antibody	100ul / \$219

21035	Anti-Cdc42/Rac/Rho Polyclonal Antibody	100ul / \$219
21036	Anti-Rab5A Rabbit Polyclonal Antibody	100ul / \$219
21037	Anti-Rab5B Rabbit Polyclonal Antibody	100ul / \$219
21038	Anti-Rab5C Rabbit Polyclonal Antibody	100ul / \$219
21039	Anti-Rab2A Rabbit Polyclonal Antibody	100ul / \$219
21041	Anti-Rab3A Rabbit Polyclonal Antibody	100ul / \$219
21042	Anti-Rab3B Rabbit Polyclonal Antibody	100ul / \$219
21043	Anti-Rab3C Rabbit Polyclonal Antibody	100ul / \$219
21044	Anti-Rab3D Rabbit Polyclonal Antibody	100ul / \$219
21046	Anti-Rab13 Rabbit Polyclonal Antibody	100ul / \$219
21047	Anti-Rab14 Rabbit Polyclonal Antibody	100ul / \$219
21049	Anti-Rab17 Rabbit Polyclonal Antibody	100ul / \$219
21050	Anti-Rab18 Rabbit Polyclonal Antibody	100ul / \$219
21054	Anti-Rab22a Rabbit Polyclonal Antibody	100ul / \$219
21056	Anti-Rab25 Rabbit Polyclonal Antibody	100ul / \$219
21058	Anti-Rab26 Rabbit Polyclonal Antibody	100ul / \$219
21061	Anti-Rab28 Rabbit Polyclonal Antibody	100ul / \$219
21062	Anti-Rab30 Rabbit Polyclonal Antibody	100ul / \$219
21063	Anti-Rab31 Rabbit Polyclonal Antibody	100ul / \$219
21068	Anti-Rab6A Rabbit Polyclonal Antibody	100ul / \$219
21069	Anti-Rab7A Rabbit Polyclonal Antibody	100ul / \$219
21071	Anti-Rab8A Rabbit Polyclonal Antibody	100ul / \$219
21073	Anti-Rab9A Rabbit Polyclonal Antibody	100ul / \$219

21076	Anti-Rab12 Rabbit Polyclonal Antibody	100ul / \$219
21078	Anti-Rab35 Rabbit Polyclonal Antibody	100ul / \$219
21081	Anti-Rab38 Rabbit Polyclonal Antibody	100ul / \$219
21084	Anti-Rab41 Rabbit Polyclonal Antibody	100ul / \$219
21090	Anti-Sar1A Rabbit Polyclonal Antibody	100ul / \$219
21097	Anti-Ran Rabbit Polyclonal Antibody	100ul / \$219
21098	Anti-RheB Rabbit Polyclonal Antibody	100ul / \$219
21110	Anti-kim-1 Rabbit Polyclonal Antibody	100ul / \$219
21112	Anti-albumin Rabbit Polyclonal Antibody	100ul / \$219
21113	Anti-ApoA1 Rabbit Polyclonal Antibody	100ul / \$219
21116	Anti-Transferrin Rabbit Polyclonal Antibody	100ul / \$219
21123	Anti-MUC1 Rabbit Polyclonal Antibody	100ul / \$219
21126	Anti-Septin9 Rabbit Polyclonal Antibody	100ul / \$219
21150	Anti-Ral Rabbit Polyclonal Antibody	100ul / \$219
21151	Anti-Rab5 Rabbit Polyclonal Antibody	100ul / \$219
21156	Anti-Rap1 Rabbit Polyclonal Antibody	100ul / \$219
21157	Anti-Rab11 Rabbit Polyclonal antibody	100ul / \$219
21158	Anti-Rab4A Rabbit polyclonal antibody	100ul / \$219
21159	Anti-Ras Rabbit polyclonal antibody	100ul / \$219
21208	Anti-STK15 Rabbit Polyclonal Antibody	100ul / \$219
21210	Anti-Cyclin B1 Rabbit Polyclonal Antibody	100ul / \$219
21218	Anti-GRB7 Rabbit Polyclonal Antibody	100ul / \$219
21220	Anti-Ras Rabbit Polyclonal Antibody	100ul / \$219

21224	Anti-BCL2 Rabbit Polyclonal Antibody	100ul / \$219
21230	Anti-Bag1 Rabbit Polyclonal Antibody	100ul / \$219
21231	Anti-Bag1 Rabbit Polyclonal Antibody	100ul / \$219
21272	Anti-EGFR Rabbit Polyclonal Antibody	100ul / \$219
21277	Anti-IDH1 Rabbit Polyclonal Antibody	100ul / \$219
21278	Anti-IDH2 Rabbit Polyclonal Antibody	100ul / \$219
21281	Anti-STK11 Rabbit Polyclonal Antibody	100ul / \$219
21283	Anti-Arf1 Rabbit Polyclonal Antibody	100ul / \$219
21013	Anti-Gα12 Rabbit Polyclonal Antibody	100ul / \$219
21006	Anti-Gα i Rabbit Polyclonal Antibody	100ul / \$219
21014	Anti-Gαq Rabbit Polyclonal Antibody	100ul / \$219
21082	Anti-Rab39 Rabbit Polyclonal Antibody	100ul / \$219
10122	Arf1 Protein Δ17 mutant	10 μg / \$347
10122-1	Arf1 Protein Δ17 mutant	50 μg / \$1386
10123-1	Arf1 Protein Δ17 Q71L mutant	50 μg / \$1386
10123	Arf1 Protein Δ17 Q71L mutant	10 μg / \$347
10124	Arf6 Protein	100 μg / \$239
10124-1	Arf6 Protein	200 μg / \$413
10125	Arf6 Protein Δ12 Q67L mutant	10 μg / \$347
10125-1	Arf6 Protein Δ12 Q67L mutant	50 μg / \$1386
10156	Arl1 Protein Q71L mutant	25 μg / \$248
10156-1	Arl1 Protein Q71L mutant	100 μg / \$605
10152	Arl3 Protein	100 μg / \$239

10152-1	Arl3 Protein	200 µg / \$413
10166	Arl3 Protein Q71L mutant	25 µg / \$248
10166-1	Arl3 Protein Q71L mutant	100 µg / \$605
10107	Cdc42 Protein	100 µg / \$239
10107-1	Cdc42 Protein	200 µg / \$413
10109	Cdc42 Protein G12V mutant	25 µg / \$248
10109-1	Cdc42 Protein G12V mutant	100 µg / \$605
10157	CST3 Protein	100 µg / \$239
10157-1	CST3 Protein	200 µg / \$413
10230	EGFR protein, L858R	100ul / \$347
10150	Fascin1 Protein	100 µg / \$239
10150-1	Fascin1 Protein	200 µg / \$413
10138	Gα13 Protein	100 µg / \$239
10138-1	Gα13 Protein	200 µg / \$413
10146	Gαi1 Protein	100 µg / \$239
10146-1	Gαi1 Protein	200 µg / \$413
10147	Gαi2 Protein Q205L mutant	10 µg / \$347
10147-1	Gαi2 Protein Q205L mutant	50 µg / \$1386
10126	Rab11 Protein	10 µg / \$347
10126-1	Rab11 Protein	50 µg / \$1386
10133	Rab35 Protein Q67L mutant	25 µg / \$248
10133-1	Rab35 Protein Q67L mutant	100 µg / \$605
10129-1	Rab3A Protein Q81L mutant	100 µg / \$605

10129	Rab3A Protein Q81L mutant	25 µg / \$248
10131	Rab4A Protein Q67L mutant	25 µg / \$248
10131-1	Rab4A Protein Q67L mutant	100 µg / \$605
10116	Rab5 Protein	100 µg / \$239
10116-1	Rab5 Protein	200 µg / \$413
10117	Rab5 Protein Q79L mutant	25 µg / \$248
10117-1	Rab5 Protein Q79L mutant	100 µg / \$605
10167	Rab7 Protein	100 µg / \$239
10167-1	Rab7 Protein	200 µg / \$413
10135	Rab7 Protein Q67L mutant	25 µg / \$248
10135-1	Rab7 Protein Q67L mutant	100 µg / \$605
10101	Rac1 Protein	100 µg / \$239
10101-1	Rac1 Protein	200 µg / \$413
10169	RalB Protein	100 µg / \$239
10169-1	RalB Protein	200 µg / \$413
10112	Ran Protein	100 µg / \$239
10112-1	Ran Protein	200 µg / \$413
10113-1	Ran Protein Q69L mutant	100 µg / \$605
10113	Ran Protein Q69L mutant	25 µg / \$248
10120	Rap1B Protein	100 µg / \$239
10120-1	Rap1B Protein	200 µg / \$413
10110	Ras Protein	100 µg / \$239
10110-1	Ras Protein	100 µg / \$413

10111	Ras Protein Q61L mutant	25 µg / \$248
10111-1	Ras Protein Q61L mutant	100 µg / \$605
10104	RhoA Protein	100 µg / \$239
10104-1	RhoA Protein	200 µg / \$413
10114	Sar1a Protein	100 µg / \$239
10114-1	Sar1a Protein	200 µg / \$413
10451	NFE2L2 Protein D29G mutant	100ul / \$347
10446	NFE2L2 Protein D77V mutant	100ul / \$347
10450	NFE2L2 Protein E79K mutant	100ul / \$347
10449	NFE2L2 Protein E79Q mutant	100ul / \$347
10448	NFE2L2 Protein L30F mutant	100ul / \$347
10447	NFE2L2 Protein R34Q mutant	100ul / \$347
10452	NFE2L2 Protein V32T mutant	100ul / \$347
10445	NFE2L2 Protein	100ul / \$347
80301	Gαi Activation Assay Kit	20 Assays / \$637
80401	Gα13 Activation Assay Kit	20 Assays / \$637
80501	Rac Activation Assay Kit	20 Assays / \$637
80601	RhoA Activation Assay Kit	20 Assays / \$637
80701	Cdc42 Activation Assay Kit	20 Assays / \$637
80801	Gαs Activation Assay Kit	20 Assays / \$637
80901	Gαo Activation Assay Kit	20 Assays / \$637
81001	Gαz Activation Assay Kit	20 Assays / \$637
81101	Ras Activation Assay Kit	20 Assays / \$637

81201	RheB Activation Assay Kit	20 Assays / \$637
81301	Rab5A Activation Assay Kit	20 Assays / \$637
81401	Rap1 Activation Assay Kit	20 Assays / \$637
81501	RalB Activation Assay Kit	20 Assays / \$637
81701	Ran Activation Assay Kit	20 Assays / \$637
81801	Sar1 Activation Assay Kit	20 Assays / \$637
82401	Arf6 Activation Assay Kit	20 Assays / \$637
82501	Rab7 Activation Assay Kit	20 Assays / \$637
82601	Rab3 Activation Assay Kit	20 Assays / \$637
82701	Rab4A Activation Assay Kit	20 Assays / \$637
82801	Rab35 Activation Assay Kit	20 Assays / \$637
82901	Arl1 Activation Assay Kit	20 Assays / \$637
83001	Arl3 Activation Assay Kit	20 Assays / \$637
83201	Rab11 Activation Assay Kit	20 Assays / \$637
83401	Rab5B Activation Assay Kit	20 Assays / \$637
83501	Rab5C Activation Assay Kit	20 Assays / \$637
83601	RalA Activation Assay Kit	20 Assays / \$637
83701	Rab5 Activation Assay Kit	20 Assays / \$637
80101	Monoclonal Anti-cGMP EIA Kit	1x96 well / \$385
80102	Monoclonal Anti-cGMP EIA Kit	5x96 well / \$1595
80103	Monoclonal Anti-cGMP EIA Kit without Acetylation	1x96 well / \$385
80104	Monoclonal Anti-cGMP EIA Kit without Acetylation	5x96 well / \$1595
80201	Monoclonal Anti-cAMP EIA Kit	1x96 well / \$385

80202	Monoclonal Anti-cAMP EIA Kit	5x96 well / \$1595
80203	Monoclonal Anti-cAMP EIA Kit without Acetylation	1x96 well / \$385
80204	Monoclonal Anti-cAMP EIA Kit without Acetylation	5x96 well / \$1595