

EPPENDORF® CONICAL TUBES

Eppendorf™ Conical Tubes with Screw Cap offer reliable sample containment and leak-free storage. SafeCode variants enhance sample identification, while the BioBased options, manufactured using renewable resources, provide an eco-friendly alternative without compromising performance.

- Ensures leak-free storage and protects sample integrity.
- Suited for a wide range of laboratory tasks, enhancing flexibility.
- High-quality materials for long-lasting, robust performance.
- Free of plasticizers, biocides or mold release agents.
- Clear, visible graduations for quick volume assessment.
- Facilitates optimal sample recovery and efficient centrifugation.
- Available in 15 mL, 25 mL, and 50 mL formats.
- 2D Datamatrix (SafeCode) option features 2-level coding for safe sample identification. ECC200 error correction function enables safe code reading with up to 25% damaged codes.

Product Number: Eppendorf™ 0030122186; Thermo Fisher 14-282-311

Closure Material: HDPE

Capacity (Metric): 50 mL

No. per Pack: 12 Racks x 25 Tubes

Quantity: 300 Pk.

Autoclavable: Not Autoclavable

Bottom Shape: Conical

Diameter (Metric): 30 mm

Purity or Quality Grade: Sterile, pyrogen-, Dnase-, Rnase-, Human and Bacterial DNA-free

Color: Clear

Material: Polypropylene

Packaging: 300 Tubes (12 Racks x 25 Tubes)

Sterility: Sterile

Clean Claims: Pyrogen, DNase, RNase, human and bacterial DNA free

Closure Type: Screw Cap

Max. RCF: 19,500 x g



Note: This product was over-purchased from Eppendorf™.