

Product Description

FISHERBRAND™ STERILE POLYSTYRENE DISPOSABLE SEROLOGICAL PIPETS WITH MAGNIFIER STRIPE - 100 ML

Fisherbrand™ Sterile Polystyrene Disposable Serological Pipets with Magnifier Stripe – 100 mL

The high-contrast black markings make the volumes easier to read, for improved optical clarity. Manufactured with optically clear polystyrene for accurate measurement. Fisherbrand™ Sterile Polystyrene Disposable Serological Pipets with Magnifier Stripe reduce reading errors and are designed for tissue culture work.

Inner bags are recyclable LDPE; outer boxes also are recyclable. Individually wrapped pipets are ideal for small-volume use requiring sterility and are available in either a clear, durable all-plastic wrap, or a color-coded paper/plastic wrap for reduced static cling.

Cat. #: Fisher Scientific 13-678-000

Calibration: To Deliver/Blow Out

Material: Polystyrene

Graduated: Yes

Negative Graduations: Yes

Tolerance: $\pm 2\%$

Color: Aqua

Type: Serological Pipet

Certifications/Compliance: Non-pyrogenic, USP XXII, USP XXII Bacterial Endotoxins Test (85); DNase/RNase Free; BSE/TSE Free

Sterility: Sterile

Graduations: 1 mL

Quantity: 10/Bag, 100/Cs.

Packaging Type: Individually Wrapped, paper-plastic wrap

Volume (Metric): 100 mL

Disposable: Single-use



Note: This product was the same Fisherbrand™ , but with 50 %discount