

BD LUER LOCK SYRINGE, 20ML

BD Luer Lock Syringe, 20mL, 48/Box 302830

- Sterile
- Single use
- Available in two tip styles: Luer-Lok™ tip and Slip tip
- DNR (Latex) free

Cat. #: BD 302830; Thermo Fisher Scientific 22-124-967

Qty./Pk: 48 / Pk

Brand: DB®

Disposable: Single-use

Graduations: 1 mL

Volume (Metric): 20 mL

Connection Type: Luer Lock

Sterility: Sterile



Note: This product was over-stocked. It is 50% off Corning list price.