

CORNING® 15 ML PP CENTRIFUGE TUBES, RACK PACKED WITH PLUG SEAL CAP, STERILE, 50/RACK, 500/CASE

Corning® 15 mL PP Centrifuge Tubes, Rack Packed with Plug Seal Cap, Sterile, 50/Rack, 500/Case

- Black printed graduations and a large white marking spot
- Maximum RCF is 12,000xg.
- Sterile, nonpyrogenic and DNase-/RNase-free
- 95 k PA (14psi) pressure tested.

Cat. #: Corning® 430052

Qty./Pk: 50 / Pk

Qty./Cs: 500 / Cs

Brand: Corning®

Capacity: 15 mL

Tube Style: Centrifuge

Bottom Style: Conical

Top Style: Threaded

Maximum RCF: 12000 xg

Graduations: Yes

Graduation Interval: 0.5 mL

Color: Clear

Cap Style: Plug seal

Cap Material: High density polyethylene

Cap Diameter: 22.9 mm (approx.)

Cap Color: Orange

Sterile: Yes

Self Standing: No

O.D.: 17.5

Length: 118.5 mm (approx.)

Feature: Disposable

Rack: Yes

Graduation Range: 1.5 – 14.5 mL



Note: This product was over-stocked from Corning. It is 50% off Corning list price.

