

## **CORNING® 15 ML PP CENTRIFUGE TUBES, RACK PACKED WITH PLUG SEAL CAP, STERILE, 50/RACK, 500/CASE**

### **Corning® 15 mL PP Centrifuge Tubes, Rack Packed with Plug Seal Cap, Sterile, 50/Rack, 500/Case**

- Black printed graduations and a large white marking spot
- Maximum RCF is 12,000xg.
- Sterile, nonpyrogenic and DNase-/RNase-free
- 95 k PA (14psi) pressure tested.

**Cat. #:** Corning® 430052

**Qty./Pk:** 50 / Pk

**Qty./Cs:** 500 / Cs

**Brand:** Corning®

**Capacity:** 15 mL

**Tube Style:** Centrifuge

**Bottom Style:** Conical

**Top Style:** Threaded

**Maximum RCF:** 12000 xg

**Graduations:** Yes

**Graduation Interval:** 0.5 mL

**Color:** Clear

**Cap Style:** Plug seal

**Cap Material:** High density polyethylene

**Cap Diameter:** 22.9 mm (approx.)

**Cap Color:** Orange

**Sterile:** Yes

**Self Standing:** No

**O.D.:** 17.5

**Length:** 118.5 mm (approx.)

**Feature:** Disposable

**Rack:** Yes

**Graduation Range:** 1.5 – 14.5 mL



**Note:** This product was over-stocked from Corning. It is 50% off Corning list price.

