

TM4SF1 MAB (AGX101 MAB)

Cat.#: 28151

Product Name: Anti-TM4SF1 (AGX101) Monoclonal Antibody

Synonyms: M3S1;TAAL6

Description: Anti-TM4SF1 (AGX101) Monoclonal Antibody

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Applications: Flow Cyt

Recommended Dilutions: Flow Cyt 1:100

Host Species: Humanized

Isotype: IgG1

Purification: Purified from cell culture supernatant by affinity chromatography

Species Reactivity: Human TM4SF1

Constituents: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with 1 µg/mL Anti-TM4SF1 (AGX101 biosimilar) mAb (28151) on Expi293 cells transfected with Human TM4SF1 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).

