

ABCE1 RABBIT MAB

Cat.#: N261609

Product Name: Anti-ABCE1 Rabbit Monoclonal Antibody

Synonyms: ABC38; ABCE1; HuHP68; OABP; RLI; RNASEL1; RNASEL2; RNS4I

UNIPROT ID: P61221

Background: Antagonizes the binding of 2-5A (5'-phosphorylated 2',5'-linked oligoadenylates) by RNase L through direct interaction with RNase L and therefore inhibits its endoribonuclease activity.

Immunogen: A synthetic peptide of human ABCE1

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R09-8C3

MW: Calculated MW: 67 kDa; Observed MW: 67 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of ABCE1 in K562, C6, 3T3, HeLa lysates using ABCE1 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
