

ABCG1 RABBIT PAB

Cat.#: N225604

Product Name: Anti-ABCG1 Rabbit pAb

Synonyms: ATP-binding cassette transporter 8; White protein homolog; ABC8; WHT1; ABCG1; ABC transporter 8; WHITE1

UNIPROT ID: P45844

Background: Functions to regulate cholesterol and phospholipid transport in macrophages. ABCG1 is highly expressed in several tissues, including brain, spleen, lung and placenta, and has been localized to the cell surface and intracellular compartments of cholesterol-laden macrophages.

Immunogen: A synthesized peptide derived from human ABCG1

Applications: WB,IHC-P,ICC/IF

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 76 kDa; Observed MW: 110 kDa

Isotype: IgG

Purification: Affinity Chromatography

Species Reactivity: Human,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Neuroscience

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of ABCG1 in HUVEC lysates using ABCG1 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
