

ACADS RABBIT MAB

Cat.#: N261785

Product Name: Anti-ACADS Rabbit Monoclonal Antibody

Synonyms: SCAD; ACAD3

UNIPROT ID: P16219

Background: Introduces a double bond at position 2 in saturated acyl-CoA's of short chain length, i.e. less than 6 carbon atoms.

Immunogen: A synthetic peptide of human ACADS

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R08-4C2

MW: Calculated MW: 44 kDa; Observed MW: 44 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of ACADS in Jurkat, C6, HeLa lysates using ACADS antibody.



Western blot analysis of ACADS in mouse liver, 3T3 lysates using ACADS antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
