

ADENOSINE A1 RECEPTOR RABBIT PAB

Cat.#: N225859

Product Name: Anti-Adenosine A1 Receptor Rabbit pAb

Synonyms: ADORA1; Adenosine receptor A1

UNIPROT ID: P30542

Background: Receptor for adenosine. The activity of this receptor is mediated by G proteins which inhibit adenylyl cyclase.

Immunogen: A synthesized peptide derived from human ADORA1

Applications: WB,IHC-P,ICC/IF,IP,FC

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20 FC: 1/50-1/100

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 37 kDa; Observed MW: 37 kDa

Isotype: IgG

Purification: Affinity Chromatography

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Neuroscience Adenosine A1

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of ADORA1 in A549 lysates using Adenosine A1 Receptor antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
