

AMBN RABBIT PAB

Cat.#: N225497

Product Name: Anti-AMBN Rabbit pAb

Synonyms: AIF; Ameloblastin

UNIPROT ID: Q9NP70

Background: Involved in the mineralization and structural organization of enamel.

Immunogen: The antiserum was produced against synthesized peptide derived from human Ameloblastin. AA range:275-324

Applications: WB,ELISA

Recommended Dilutions: WB: 1/500-1/1000 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 48 kDa; Observed MW: 48 kDa

Isotype: IgG

Purification: Affinity Chromatography

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of AMBN in HeLa lysates using Ameloblastin antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
