

ANKRD20A1 RABBIT PAB

Cat.#: N225201

Product Name: Anti-ANKRD20A1 Rabbit pAb

Synonyms: ANKRD20A1; ANKRD20A; Ankyrin repeat domain-containing protein 20A1

UNIPROT ID: Q5TYW2

Background: Located in plasma membrane.

Immunogen: Synthesized peptide derived from the Internal region of human A20.

Applications: WB,ELISA

Recommended Dilutions: WB: 1/500-1/1000 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 94 kDa; Observed MW: 94 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of ANKRD20A1 in 3T3 lysates using ANKRD20A1 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
