

ARFRP1 RABBIT MAB

Cat.#: N261581

Product Name: Anti-ARFRP1 Rabbit Monoclonal Antibody

Synonyms: ARP; Arp1; ARL18

UNIPROT ID: Q13795

Background: Trans-Golgi-associated GTPase that regulates protein sorting.

Immunogen: A synthetic peptide of human ARFRP1

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R08-5C5

MW: Calculated MW: 23 kDa; Observed MW: 23 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of ARFRP1/ARP in C6 lysates using ARFRP1 antibody.