

## ARHGAP25 RABBIT MAB

**Cat.#:** N261869

**Product Name:** Anti-ARHGAP25 Rabbit Monoclonal Antibody

**Synonyms:** Rho GTPase activating protein 25; KAIA0053

**UNIPROT ID:** P42331

**Background:** GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state.

**Immunogen:** Recombinant protein of human ARHGAP25

**Applications:** WB,IHC-P

**Recommended Dilutions:** WB: 1/500-1/1000 IHC: 1/50-1/100

**Host Species:** Rabbit

**Clonality:** Rabbit Monoclonal

**Clone ID:** R07-6G6

**MW:** Calculated MW: 73 kDa; Observed MW: 73 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human,Rat


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Signal Transduction

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of ARHGAP25 in rat brain, Raoms, Jurkat lysates using ARHGAP25 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---