

AVEN RABBIT MAB

Cat.#: N261899

Product Name: Anti-AVEN Rabbit Monoclonal Antibody

Synonyms: PDCD12

UNIPROT ID: Q9NQS1

Background: Protects against apoptosis mediated by Apaf-1.

Immunogen: A synthetic peptide of human AVEN

Applications: WB,IHC-P

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R06-7A9

MW: Calculated MW: 39 kDa; Observed MW: 50 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of AVEN in rat Brain, MCF-7, 3T3, HeLa lysates using AVEN antibody.



Immunohistochemistry analysis of paraffin-embedded mouse testis using AVEN antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.