

BAG2 RABBIT MAB

Cat.#: N261576

Product Name: Anti-BAG2 Rabbit Monoclonal Antibody

Synonyms: BAG 2; Bag2

UNIPROT ID: O95816

Background: Inhibits the chaperone activity of HSP70/HSC70 by promoting substrate release.

Immunogen: A synthetic peptide of human BAG2

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R09-7C3

MW: Calculated MW: 24 kDa; Observed MW: 24 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of BAG2 in Jurkat, HeLa lysates using BAG2 antibody.