

BRPF3 RABBIT PAB

Cat.#: S218374

Product Name: Anti-BRPF3 Rabbit Polyclonal Antibody

Synonyms:

UNIPROT ID: Q9ULD4 (Gene Accession - BC117387)

Background: Component of the MOZ/MORF complex which has a histone H3 acetyltransferase activity.

Immunogen: Fusion protein of human BRPF3

Applications: ELISA, IHC

Recommended Dilutions: IHC: 25-100; ELISA: 5000-10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Epigenetics and Nuclear Signaling, Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin embedded Human brain tissue using 218374(BRPF3 Antibody) at a dilution of 1/40(Cytoplasm and Nucleus).



In comparison with the IHC on the left, the same paraffin-embedded Human brain tissue is first treated with the fusion protein and then with 218374(Anti-BRPF3 Antibody) at dilution 1/40.