

C11ORF49 RABBIT PAB

Cat.#: S212331

Product Name: Anti-C11orf49 Rabbit Polyclonal Antibody

Synonyms:

UNIPROT ID: Q9H6J7 (Gene Accession - BC009455)

Background: The function of this protein remains unknown.

Immunogen: Fusion protein of human C11orf49

Applications: ELISA, IHC

Recommended Dilutions: IHC: 25-100; ELISA: 5000-10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human, Mouse, Rat


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin embedded Human thyroid cancer tissue using 212331(C11orf49 Antibody) at a dilution of 1/20(Cytoplasm).



In comparison with the IHC on the left, the same paraffin-embedded Human thyroid cancer tissue is first treated with the fusion protein and then with 212331(Anti-C11orf49 Antibody) at dilution 1/20.