

C4 BINDING PROTEIN RABBIT MAB

Cat.#: N261944

Product Name: Anti-C4 Binding Protein Rabbit Monoclonal Antibody

Synonyms: C4BP; C4BPB; PRP

UNIPROT ID: P20851

Background: Controls the classical pathway of complement activation. It binds as a cofactor to C3b/C4b inactivator (C3bINA), which then hydrolyzes the complement fragment C4b.

Immunogen: Recombinant protein of human C4 binding protein

Applications: WB,IHC-F,IHC-P,ICC/IF,IP

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R02-8F3

MW: Calculated MW: 28 kDa; Observed MW: 40 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Immunology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of C4 binding protein in Human serum lysates using C4 Binding Protein antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
