

C8 ALPHA RABBIT PAB

Cat.#: N225236

Product Name: Anti-C8 alpha Rabbit pAb

Synonyms: Complement component C8 alpha chain (Complement component 8 subunit alpha)

UNIPROT ID: P07357

Background: Constituent of the membrane attack complex (MAC) that plays a key role in the innate and adaptive immune response by forming pores in the plasma membrane of target cells.

Immunogen: Synthesized peptide derived from human C8 α . at AA range: 201-250

Applications: WB,ELISA

Recommended Dilutions: WB: 1/500-1/1000 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 65 kDa; Observed MW: 60 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg^{2+} and Ca^{2+}), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Immunology

Storage & Shipping: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing



Western blot analysis of C8 alpha in mouse kidney, mouse heart lysates using C8 alpha antibody.