

## CACTIN RABBIT PAB

**Cat.#:** N225213

**Product Name:** Anti-Cactin Rabbit pAb

**Synonyms:** CACTIN C19orf29

**UNIPROT ID:** Q8WUQ7

**Background:** Involved in the regulation of innate immune response.

**Immunogen:** The antiserum was produced against synthesized peptide derived from the Internal region of human CACTIN. AA range:141-190

**Applications:** WB,ELISA

**Recommended Dilutions:** WB: 1/500-1/1000 ELISA: 1/10000

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Clone ID:** -

**MW:** Calculated MW: 89 kDa; Observed MW: 95 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human,Mouse,Rat

**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Epigenetics and Nuclear Signaling

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Cactin in hepg2 lysates using Cactin antibody.