

CALRETININ RABBIT PAB

Cat.#: N225840

Product Name: Anti-Calretinin Rabbit pAb

Synonyms: Calretinin; CR; 29 kDa calbindin

UNIPROT ID: P22676

Background: Calretinin is a calcium-binding protein which is abundant in auditory neurons.

Immunogen: A synthesized peptide derived from human Calretinin

Applications: WB,IHC-P,ICC/IF,FC

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 FC: 1/50-1/100

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 32 kDa; Observed MW: 29 kDa

Isotype: IgG

Purification: Affinity Chromatography

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Tags & Cell Markers

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Calretinin in SH-SY5Y lysates using Calretinin antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
