

CAV GAMMA 5 RABBIT PAB

Cat.#: N225536

Product Name: Anti-Cav gamma 5 Rabbit pAb

Synonyms: acng5 antibody; TARP gamma-5 antibody; Transmembrane AMPAR regulatory protein gamma-5 antibody

UNIPROT ID: Q9UF02

Background: Regulates the gating properties of AMPA-selective glutamate receptors (AMPA-Rs).

Immunogen: Synthetic Peptide of L-type Ca⁺⁺ CP γ 5

Applications: WB,IHC-P

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 31 kDa; Observed MW: 35 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Cav gamma 5 in Human Brain, mouse Brain, rat Brain lysates using Cav Cav gamma 5 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
