

CD1C RABBIT PAB

Cat.#: N225175

Product Name: Anti-CD1C Rabbit pAb

Synonyms: CD1C; T-cell surface glycoprotein CD1c; CD1c

UNIPROT ID: P29017

Background: Antigen-presenting protein that binds self and non-self lipid and glycolipid antigens and presents them to T-cell receptors on natural killer T-cells.

Immunogen: The antiserum was produced against synthesized peptide derived from the Internal region of human CD1C. AA range:211-260

Applications: WB,ELISA

Recommended Dilutions: WB: 1/500-1/1000 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 38 kDa; Observed MW: 38 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human

Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Immunology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of CD1C in K562 lysates using CD1C antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
