

CD3E(MOSUNETUZUMAB) MAB

Cat.#: 28575

Product Name: Anti-CD3E(Mosunetuzumab) Monoclonal Antibody

Synonyms: CD3e, T3E

Description: Anti-CD3E(mosunetuzumab) Monoclonal Antibody

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Applications: Flow Cyt

Recommended Dilutions: Flow Cyt 1:100

Host Species: Humanized

Isotype: IgG1

Purification: Purified from cell culture supernatant by affinity chromatography

Species Reactivity: Human CD3E

Constituents: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. CD3E protein is highly expressed inside Jurkat cells. Flow cytometry analysis with 1 μ g/mL Anti-CD3E(mosunetuzumab biosimilar) mAb (28575) (Blue histogram) or isotype control mAb (Red histogram) on Jurkat cells.