

## CD3E(MOSUNETUZUMAB) MAB

**Cat.#:** 28575

**Product Name:** Anti-CD3E(Mosunetuzumab) Monoclonal Antibody

**Synonyms:** CD3e, T3E

**Description:** Anti-CD3E(mosunetuzumab ) Monoclonal Antibody

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Applications:** Flow Cyt

**Recommended Dilutions:** Flow Cyt 1:100

**Host Species:** Humanized

**Isotype:** IgG1

**Purification:** Purified from cell culture supernatant by affinity chromatography

**Species Reactivity:** Human CD3E

**Constituents:** Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. CD3E protein is highly expressed inside Jurkat cells. Flow cytometry analysis with 1 µg/mL Anti-CD3E(mosunetuzumab biosimilar) mAb (28575) (Blue histogram) or isotype control mAb (Red histogram) on Jurkat cells.