

## CD3G RABBIT MAB

**Cat.#:** N262001

**Product Name:** Anti-CD3G Rabbit Monoclonal Antibody

**Synonyms:** CD3 gamma; CD3g; T3G

**UNIPROT ID:** P09693

**Background:** The CD3 complex mediates signal transduction.

**Immunogen:** A synthetic peptide of human CD3G

**Applications:** WB, ICC/IF

**Recommended Dilutions:** WB: 1/500-1/1000 IF: 1/50-1/200

**Host Species:** Rabbit

**Clonality:** Rabbit Monoclonal

**Clone ID:** R02-4D1

**MW:** Calculated MW: 20 kDa; Observed MW: 20 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Immunology

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of CD3G in Jurkat lysates using CD3G antibody.



Immunocytochemistry analysis of CD3G (green) in Jurkat using CD3G antibody, and DAPI (blue).



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---