

## CD99 RABBIT MAB

**Cat.#:** N262008

**Product Name:** Anti-CD99 Rabbit Monoclonal Antibody

**Synonyms:** MIC2; HBA71; MIC2X; MIC2Y; MSK5X

**UNIPROT ID:** P14209

**Background:** Involved in T-cell adhesion processes. It is involved in spontaneous rosette formation with erythrocytes.

**Immunogen:** A synthetic peptide of human CD99

**Applications:** WB,IHC-P,IP

**Recommended Dilutions:** WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

**Host Species:** Rabbit

**Clonality:** Rabbit Monoclonal

**Clone ID:** R06-8G3

**MW:** Calculated MW: 19 kDa; Observed MW: 30 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human,Mouse


**Conjugation:** Unconjugated

**Modification:** Unmodified


**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Immunology


**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using CD99 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of CD99 in MOLT4 lysates using CD99 antibody.



Western blot analysis of CD99 in 3T3 lysates using CD99 antibody.