

CDKN2AIP RABBIT PAB

Cat.#: S216699

Product Name: Anti-CDKN2AIP Rabbit Polyclonal Antibody

Synonyms: CARF

UNIPROT ID: Q9NXV6 (Gene Accession - BC022270)

Background: CDKN2AIP (CDKN2A-interacting protein), also known as CARF, is a 580 amino acid protein that activates p53 via p14 ARF (alternate reading frame)-dependent and independent pathways. CDKN2AIP-dependent activation of p53, a protein that up-regulates growth arrest and apoptosis-related genes in response to stress signals, leads to an enhancement of p53 function. Expression levels of CDKN2AIP and p53 show an inverse relationship that is caused by a negative-feedback control via a proteasome-mediated degradation pathway. CDKN2AIP is expressed ubiquitously across tissue samples and, along with p14 ARF, is localized to the perinucleolar region within the nucleus.

Immunogen: Fusion protein of human CDKN2AIP

Applications: ELISA, IHC

Recommended Dilutions: IHC: 25-100; ELISA: 2000-5000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Epigenetics and Nuclear Signaling, Cancer

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin embedded Human brain tissue using 216699(CDKN2AIP Antibody) at a dilution of 1/20(Cytoplasm).



In comparison with the IHC on the left, the same paraffin-embedded Human brain tissue is first treated with the fusion protein and then with 216699(Anti-CDKN2AIP Antibody) at dilution 1/20.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
