

CLEAVED-CASPASE 9 RABBIT PAB

Cat.#: N225350

Product Name: Anti-Cleaved-Caspase 9 Rabbit pAb

Synonyms: CASP9; MCH6; Caspase-9; CASP-9; Apoptotic protease Mch-6; Apoptotic protease-activating factor 3; APAF-3; ICE-like apoptotic protease 6; ICE-LAP6

UNIPROT ID: Q9JHK1

Background: Involved in the activation cascade of caspases responsible for apoptosis execution. Binding of caspase-9 to Apaf-1 leads to activation of the protease which then cleaves and activates caspase-3.

Immunogen: The antiserum was produced against synthesized peptide derived from rat Caspase 9. AA range:323-372

Applications: ICC/IF,WB,IHC-F,IHC-P,ELISA

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200
ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 50,17 kDa; Observed MW: 17 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Mouse,Rat


Conjugation: Unconjugated

Modification: Cleaved

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Cleaved-Caspase 9 in various lysates using Cleaved-Caspase 9 antibody.