

CLPTM1 RABBIT MAB

Cat.#: N262047

Product Name: Anti-CLPTM1 Rabbit Monoclonal Antibody

Synonyms: Cleft lip and palate transmembrane protein 1

UNIPROT ID: O96005

Background: May play a role in T-cell development.

Immunogen: A synthetic peptide of human CLPTM1

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R05-316

MW: Calculated MW: 76 kDa; Observed MW: 90 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of CLPTM1 in 293T, U87-MG, Jurkat lysates using CLPTM1 antibody.