

COILIN RABBIT PAB

Cat.#: N225985

Product Name: Anti-Coilin Rabbit pAb

Synonyms: CLN80; COIL; Coilin; Coilin p80

UNIPROT ID: P38432

Background: Is a component of the nuclear coiled bodies (CBS) which are involved in the function or assembly/disassembly of nucleoplasmic snRNPs. During mitosis, CBS disassemble, coinciding with a mitotic-specific phosphorylation of p80 coilin.

Immunogen: A synthesized peptide derived from human Coilin

Applications: WB,FC

Recommended Dilutions: WB: 1/500-1/1000 FC: 1/50-1/100

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 63 kDa; Observed MW: 80 kDa

Isotype: IgG

Purification: Affinity Chromatography

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Neuroscience

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Coilin in Jurkat lysates using Coilin antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
