

## COLLAGEN VI ALPHA 1 RABBIT MAB

**Cat.#:** N262053

**Product Name:** Anti-Collagen VI alpha 1 Rabbit Monoclonal Antibody

**Synonyms:** OPLL; BTHLM1; UCHMD1

**UNIPROT ID:** P12109

**Background:** Collagen VI acts as a cell-binding protein.

**Immunogen:** Recombinant protein of human Collagen VI

**Applications:** WB,IHC-F,IHC-P,ICC/IF

**Recommended Dilutions:** WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

**Host Species:** Rabbit

**Clonality:** Rabbit Monoclonal

**Clone ID:** R02-819

**MW:** Calculated MW: 109 kDa; Observed MW: 147 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Cell Biology

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Collagen VI in Human placenta lysates using Collagen VI alpha 1 antibody.