

## COLLAGEN XIX ALPHA 1 RABBIT PAB

**Cat.#:** N225231

**Product Name:** Anti-Collagen XIX alpha 1 Rabbit pAb

**Synonyms:** Collagen alpha-1(XIX) chain (Collagen alpha-1(Y) chain)

**UNIPROT ID:** Q14993

**Background:** May act as a cross-bridge between fibrils and other extracellular matrix molecules.

**Immunogen:** Synthesized peptide derived from human COL19A1. at AA range: 161-210

**Applications:** WB,ELISA

**Recommended Dilutions:** WB: 1/500-1/1000 ELISA: 1/10000

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Clone ID:** -

**MW:** Calculated MW: 115 kDa; Observed MW: 130 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Cell Biology

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Collagen XIX alpha 1 in HELA lysates using Collagen XIX alpha 1 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---