

COLLAGEN XXVIII ALPHA 1 RABBIT PAB

Cat.#: N225209

Product Name: Anti-Collagen XXVIII alpha 1 Rabbit pAb

Synonyms: Collagen alpha-1(XXVIII) chain

UNIPROT ID: Q2UY09

Background: May act as a cell-binding protein.

Immunogen: The antiserum was produced against synthesized peptide derived from the N-terminal region of human COL28A1. AA range:111-160

Applications: WB,ELISA

Recommended Dilutions: WB: 1/500-1/1000 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 117 kDa; Observed MW: 117 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse

Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Collagen XXVIII alpha 1 in 3T3 HEPG2 lysates using Collagen XXVIII alpha 1 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
