

CRIPTO 3 RABBIT PAB

Cat.#: N225205

Product Name: Anti-Cripto 3 Rabbit pAb

Synonyms: TDGFIP3; CRIPTO3; TDGF2; TDGF3; Putative teratocarcinoma-derived growth factor 3; Cripto-3 growth factor; Epidermal growth factor-like cripto protein CR3; Teratocarcinoma-derived growth factor 1 pseudogene 3

UNIPROT ID: P51864

Background: Could play a role in the determination of the epiblastic cells that subsequently give rise to the mesoderm.

Immunogen: The antiserum was produced against synthesized peptide derived from the Internal region of human TDGFIP3. AA range:1-50

Applications: WB,ELISA

Recommended Dilutions: WB: 1/500-1/1000 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 21 kDa; Observed MW: 21 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Cripto 3 in HepG2, 293T , AD293 , HeLa lysates using Cripto 3 antibody.