

DAPI1 RABBIT MAB

Cat.#: N262112

Product Name: Anti-DAPI1 Rabbit Monoclonal Antibody

Synonyms: Death-associated protein 1; DAP-1

UNIPROT ID: P51397

Background: Negative regulator of autophagy. Involved in mediating interferon-gamma-induced cell death.

Immunogen: A synthetic peptide of human DAPI1

Applications: WB,IHC-F,IHC-P,ICC/IF,IP

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R05-2J5

MW: Calculated MW: 11 kDa; Observed MW: 11 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of DAPI1 in MCF-7, C6, 3T3 lysates using DAPI1 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
