

DESMOCOLLIN 3 RABBIT MAB

Cat.#: N261610

Product Name: Anti-Desmocollin 3 Rabbit Monoclonal Antibody

Synonyms: Cadherin family member 3; Desmocollin-4; HT-CP

UNIPROT ID: Q14574

Background: Component of intercellular desmosome junctions.

Immunogen: A synthetic peptide of human Desmocollin 3

Applications: WB, ICC/IF

Recommended Dilutions: WB: 1/500-1/1000 IF: 1/50-1/200

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R06-4G2

MW: Calculated MW: 100 kDa; Observed MW: 100 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Desmocollin 3 in K562, HeLa lysates using Desmocollin 3 antibody.



Immunocytochemistry analysis of Desmocollin 3 (green) in HeLa using Desmocollin 3 antibody, and DAPI (blue).



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
