

DHRS7B RABBIT PAB

Cat.#: S213023

Product Name: Anti-DHRS7B Rabbit Polyclonal Antibody

Synonyms: CGI-93; SDR32C1

UNIPROT ID: Q6IAN0 (Gene Accession - BC004126)

Background: This gene is located within the Smith-Magenis syndrome region on chromosome 17. It encodes a protein of unknown function.

Immunogen: Fusion protein of human DHRS7B

Applications: ELISA, IHC

Recommended Dilutions: IHC: 50-100; ELISA: 5000-10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Metabolism

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin embedded Human liver cancer tissue using 213023(DHRS7B Antibody) at a dilution of 1/80(Cytoplasm).



In comparison with the IHC on the left, the same paraffin-embedded Human liver cancer tissue is first treated with the fusion protein and then with 213023(Anti-DHRS7B Antibody) at dilution 1/80.