

DHX8 RABBIT MAB

Cat.#: N262133

Product Name: Anti-DHX8 Rabbit Monoclonal Antibody

Synonyms: DDX8; Dhr2; HRH1; PRP22; PRPF22

UNIPROT ID: Q14562

Background: Facilitates nuclear export of spliced mRNA by releasing the RNA from the spliceosome.

Immunogen: A synthetic peptide of human DHX8

Applications: WB, ICC/IF

Recommended Dilutions: WB: 1/500-1/1000 IF: 1/50-1/200

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R04-117

MW: Calculated MW: 139 kDa; Observed MW: 139 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of DHX8 in HeLa, 293T lysates using DHX8 antibody.