

DKK3 RABBIT MAB

Cat.#: N262138

Product Name: Anti-DKK3 Rabbit Monoclonal Antibody

Synonyms: DKK3; REIC; Dickkopf-related protein 3; Dickkopf-3; Dkk-3; hDkk-3

UNIPROT ID: Q9UBP4

Background: Antagonizes canonical Wnt signaling by inhibiting LRP5/6 interaction with Wnt and by forming a ternary complex with the transmembrane protein KREMEN that promotes internalization of LRP5/6.

Immunogen: Recombinant protein of human Dkk3

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R03-8G7

MW: Calculated MW: 38 kDa; Observed MW: 70 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Dkk3 in rat Brain lysates using DKK3 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
