

EMSY RABBIT PAB

Cat.#: N226026

Product Name: Anti-EMSY Rabbit pAb

Synonyms: BRCA2-interacting transcriptional repressor EMSY

UNIPROT ID: Q7Z589

Background: Regulator which is able to repress transcription, possibly via its interaction with a multiprotein chromatin remodeling complex that modifies the chromatin.

Immunogen: A synthetic peptide of human EMSY

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 141 kDa; Observed MW: 141 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of EMSY in HL-60, U2OS, C6 lysates using EMSY antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
