

## EMSY RABBIT PAB

**Cat.#:** N226026

**Product Name:** Anti-EMSY Rabbit pAb

**Synonyms:** BRCA2-interacting transcriptional repressor EMSY

**UNIPROT ID:** Q7Z589

**Background:** Regulator which is able to repress transcription, possibly via its interaction with a multiprotein chromatin remodeling complex that modifies the chromatin.

**Immunogen:** A synthetic peptide of human EMSY

**Applications:** WB,IP

**Recommended Dilutions:** WB: 1/500-1/1000 IP: 1/20

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Clone ID:** -

**MW:** Calculated MW: 141 kDa; Observed MW: 141 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human,Rat


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Cell Biology

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of EMSY in HL-60, U2OS, C6 lysates using EMSY antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---