


**EYA2 RABBIT MAB****Cat.#:** N263323**Product Name:** Anti-EYA2 Rabbit Monoclonal Antibody**Synonyms:** EAB1**UNIPROT ID:** O00167**Background:** Functions both as protein phosphatase and as transcriptional coactivator for SIX1, and probably also for SIX2, SIX4 and SIX5.**Immunogen:** A synthetic peptide of human EYA2**Applications:** WB,IP**Recommended Dilutions:** WB: 1/500-1/1000 IP: 1/20**Host Species:** Rabbit**Clonality:** Rabbit Monoclonal**Clone ID:** R01-6A1**MW:** Calculated MW: 59 kDa; Observed MW: 59 kDa**Isotype:** IgG**Purification:** Affinity Purified**Species Reactivity:** Human**Conjugation:** Unconjugated**Modification:** Unmodified**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide**Research Areas:** Epigenetics and Nuclear Signaling**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing

Western blot analysis of EYA2 in Hela, A549, U2OS lysates using EYA2 antibody.

