

EYA2 RABBIT MAB

Cat.#: N263323

Product Name: Anti-EYA2 Rabbit Monoclonal Antibody

Synonyms: EAB1

UNIPROT ID: O00167

Background: Functions both as protein phosphatase and as transcriptional coactivator for SIX1, and probably also for SIX2, SIX4 and SIX5.

Immunogen: A synthetic peptide of human EYA2

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R01-6A1

MW: Calculated MW: 59 kDa; Observed MW: 59 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of EYA2 in HeLa, A549, U2OS lysates using EYA2 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
