


FAS LIGAND RABBIT PAB**Cat.#:** N225361**Product Name:** Anti-Fas Ligand Rabbit pAb**Synonyms:** FASLG; APT1LG1; CD95L; FASL; TNFSF6; Tumor necrosis factor ligand superfamily member 6; Apoptosis antigen ligand; APTL; CD95 ligand; CD95-L; Fas antigen ligand; Fas ligand; FasL; CD antigen CD178**UNIPROT ID:** P48023**Background:** Cytokine that binds to TNFRSF6/FAS, a receptor that transduces the apoptotic signal into cells.**Immunogen:** The antiserum was produced against synthesized peptide derived from human FAS ligand. AA range:101-150**Applications:** WB,IHC-F,IHC-P,ICC/IF,ELISA**Recommended Dilutions:** WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200
ELISA: 1/10000**Host Species:** Rabbit**Clonality:** Rabbit Polyclonal**Clone ID:** -**MW:** Calculated MW: 31 kDa; Observed MW: 31 kDa**Isotype:** IgG**Purification:** Affinity Purified**Species Reactivity:** Human,Mouse**Conjugation:** Unconjugated**Modification:** Unmodified**Constituents:** PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide**Research Areas:** Cell Biology**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Fas Ligand in various lysates using Fas Ligand antibody.