


FETUIN A RABBIT MAB**Cat.#:** N262216**Product Name:** Anti-Fetuin A Rabbit Monoclonal Antibody**Synonyms:** 16G2; A2HS; AHS; Ahsg; FETUA; FETUB; FetuinA; FetuinB; Gugu; HSGA; IRL685; PRO2743**UNIPROT ID:** P02765**Background:** Promotes endocytosis, possesses opsonic properties and influences the mineral phase of bone. Shows affinity for calcium and barium ions.**Immunogen:** A synthetic peptide of human Fetuin A**Applications:** WB,IP**Recommended Dilutions:** WB: 1/500-1/1000 IP: 1/20**Host Species:** Rabbit**Clonality:** Rabbit Monoclonal**Clone ID:** R03-5A4**MW:** Calculated MW: 39 kDa; Observed MW: 45-50 kDa**Isotype:** IgG**Purification:** Affinity Purified**Species Reactivity:** Human,Mouse**Conjugation:** Unconjugated**Modification:** Unmodified**Constituents:** PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide**Research Areas:** Cardiovascular**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing

Western blot analysis of Fetuin A in HepG2, 3T3, U937 lysates using Fetuin A antibody.

