

## FGF11 RABBIT PAB

**Cat.#:** N225191

**Product Name:** Anti-FGF11 Rabbit pAb

**Synonyms:** FGF11; FHF3; Fibroblast growth factor 11; FGF-11; Fibroblast growth factor homologous factor 3; FHF-3

**UNIPROT ID:** Q92914

**Background:** Probably involved in nervous system development and function.

**Immunogen:** The antiserum was produced against synthesized peptide derived from the Internal region of human FGF11. AA range:21-70

**Applications:** WB,ELISA

**Recommended Dilutions:** WB: 1/500-1/1000 ELISA: 1/10000

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Clone ID:** -

**MW:** Calculated MW: 25 kDa; Observed MW: 25 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human,Mouse,Rat


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Cardiovascular

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of FGF11 in NIH3T3, HeLa lysates using FGF11 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---